

To ensure you are using the proper size boots for your horse please take a measurement of your horse's fetlock.

We recommend using a soft tape measure. If you do not have a soft tape measure you can use something like string, wrap it around the fetlock, mark where the end meets the rest of the string and then simply lay it flat next to a tape measure, ruler or yard stick to get your measurement.



Small - 8.5 - 9.5 in.
Medium - 9.5 - 10.5 in.
Large - 10.5 - 11.5 in.
X-Large - 11.5 - 12.5 in.
XX-Large- 12.5 - 13.5 in.

Our Measuring chart is good for both our Iconoclast® Orthopedic Sport Boots as well as our Iconoclast® Orthopedic/Rehabilitation Boots.